

Welcome to the Bio-Strategy Source Flyer

The best deals in Instrumentation, Consumables, Chemicals, Assays and Kits

The October to December edition features the following brands and many more.



Chemicals on sale at Bio-Strategy

SAVE \$100 ON CHEMICALS

Buy \$600 worth of mixed Laboratory Chemicals including VWR Chemicals and get \$100 off!

Bio-Strategy - your One Stop Shop for Laboratory Chemicals!



DigiDoc-It® Basic Gel Imaging System

30% OFF

DigiDoc-It Gel Imaging System with UV Transilluminator

Fluorescence and colourmetric applications with 17.9 MP digital colour camera

analytikjena
An Endress+Hauser Company



Pyrex Premium Volumetric Flasks!

40% OFF

PYREX® Volumetric Flasks, borosilicate glass with polyethylene stopper

SciLabware
AZLON® MBL® PYREX® QUICKFIT® WHEATON®



Panasonic Dual Heat Sterilisation Incubators

Panasonic now introduces Dual Heat sterilisation into their Incubator line up

Complimenting the existing innovative H₂O₂ decontamination models, Panasonic gives customers the choice for their preferred incubator decontamination method.

NEW



Panasonic

GE Healthcare Protein Xtra

UP TO 48% DISCOUNT!

NEW Protein Xtra Brochure!
Simplify your protein workflow and profit from up to 48% discount!



Distributor
GE Healthcare



CRYOEXTRACT CXT 353

CALL US FOR MORE DETAILS!

Frozen Sample Aliquoting System

Providing excellent sample integrity by preventing the freeze / thaw process.



CRYOEXTRACT

VWR Transfer Pipettes

15% OFF

Great savings on VWR Transfer Pipettes Made from LDPE, disposable

VWR
We Enable Science



Bio-Strategy PTY LTD.

T: 1300 727 696
W: bio-strategy.com
E: sales.au@bio-strategy.com



bio-strategy
delivering technology

SOURCE

GE HEALTHCARE PROTEIN XTRA

NEW Protein Xtra brochure!

Simplify your protein workflows and profit from up to 48% discount!

Select the right tools commonly used in Protein Workflows and profit from huge savings offered with Protein Xtra! Download or order your Protein Xtra brochure today.

- **Isolation and Cloning:** Including great savings on Ficoll Paque, Percoll, Illustra PCR Purification and Amplification kits and many more!
- **Sample Preparation:** Whatman Syringe filters, filter paper, membranes, Benchkote, PD-10 desalting columns
- **Protein Purification:** GST and Histidine tagged chromatography columns, Size Exclusion columns, Desalting columns, Vivaspin concentrators, filtration, Akta
- **Analysis:** Amersham ECL Western blotting, Molecular Weight Markers, Blotting Paper, Syringe filters

Great savings on GE Healthcare and Whatman products!



UP TO 48% OFF



Distributor
GE Healthcare

Ref. Code: **GE-XTRA**

Valid until 31/03/2018

Description	Unit of Measure	Code	Minimum ordering quantity to get the Special Price	Total List Price for bundle (excl. GST)	Total Special Price for Bundle (excl. GST)
Percoll 1L	1 * 1 l	GEHE17-0891-01	4	\$2,044.00	\$1,737.40
Ficoll Paque Plus 6X100ML	6 * 100 ml	GEHE17-1440-02	4	\$1,140.00	\$969.00
Desalting Column PD-10	1 * 30 items	GEHE17-0851-01	2	\$1,012.00	\$809.60
ECL Prime Western Blotting Detection Reagent for 1000cm2 membrane area	1 kit	GEHERPN2232	2	\$632.00	\$505.60
Isolation Kit - RNASPIN MINI	1 * 250 tests	GEHE25-0500-72	1	\$1,194.00	\$895.50
Amersham Full Range Rainbow ECL Molecular Weight Marker	1 * 250 µl	GEHERPN800E	1	\$315.00	\$252.00
Affinity chromatography column HiTrap™ TALON crude 5X5ML	1 * 5 items	GEHE28-9537-67	2	\$1,812.00	\$1,540.20
Syringe filters, GD/X™, Whatman™, PES,25MM 0.45UM	1 * 150 items	WHAT6876-2504	1	\$520.00	\$373.97

And many more!

EPPENDORF 5424 AND 5424R CENTRIFUGE



SAVE NOW

Equipped for all modern life science applications

Eppendorf 5424 and 5424R (refrigerated) centrifuges are the new laboratory standard. With their 24-place capacity and speed up to 21,130xg they are perfectly equipped for all life science applications in Eppendorf and PCR tubes

- 5424R centrifuge cools down to 4°C in only 8 minutes
- 5424R centrifuge temperature range from -10°C to +40°C
- Soft-touch one-finger lid closure for ergonomic operation
- OptiBowl® design for quiet operation, even without a rotor lid
- SOFT brake to protect sensitive samples
- Certified aerosol-tight rotor comes standard for safe centrifugation of hazardous samples



Ref Code: **EPPE-5424**

Valid until 31/12/2017

Description	Code	Sale Price (excl. GST)
Centrifuge 5424R refrigerated knob version with FA-45-24-11 rotor	EPPE5404000.480	\$5,945.00
Centrifuge 5424 knob version with FA-45-24-11 rotor	EPPE5424000.487	\$2,819.00

EPPENDORF ADVANTAGE PROMOTION

Save now on Eppendorf 5430R Centrifuge, Eppendorf ThermoMixer C and Eppendorf 5810R Centrifuge - Get the Eppendorf Advantage!

Eppendorf 5430R refrigerated microcentrifuge with multipurpose capabilities. Combines the best features of a microcentrifuge (small footprint) and multipurpose centrifuge (versatility) in one instrument.

The Eppendorf ThermoMixer C is for heating, cooling, and mixing almost any of your lab vessels. Reliable prevention of condensate formation in small vessels can be achieved with the heated lid Eppendorf ThermoTop (optional).

Eppendorf 5810R refrigerated centrifuge is a workhorse for medium to high-throughput laboratories. It combines extraordinary versatility and capacity for both tubes and plates with an extraordinary compact footprint.

NEW



advantage
an Eppendorf Promotion

Special Prices

Ref Code: EPPE-ADVAN

Valid until 31/12/2017

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
5430R Bundle (knob) with FA-45-30-11 and F-35-6-30 rotors and a box each of 15 and 50ml conical tubes	EPPE5428000.384	\$11,600.00	\$9,860.00
5430R Bundle (foil) with FA-45-30-11 and F-35-6-30 rotors and a box each of 15 and 50ml conical tubes	EPPE5428000.880	\$11,600.00	\$9,860.00
ThermoMixer C Tubes Bundle includes: ThermoMixer C with SmartBlocks for 1.5ml and 2.0ml tubes, 1 box each of 1.5ml and 2.0ml SafeLock tubes	EPPE2232000.083	\$5,235.80	\$4,607.50
ThermoMixer C Plates Bundle includes: ThermoMixer C with SmartBlocks for 1.5ml and Plates. 1 box of 1.5mL SafeLock tubes and a ThermoTop	EPPE2232000.084	\$5,972.80	\$5,256.06
Centrifuge 5810R, refrigerated, with Rotor S-4-104 incl. adapters for 15/50ml conical tubes, 15/50 conical tubes and two plate buckets	EPPE5811000.983	NEW!	\$16,663.00
Centrifuge 5810R, refrigerated, with Rotor S-4-104 incl. adapters for 15/50ml conical tubes	EPPE5811000.681	\$17,540.00	\$15,786.00

MOLECULAR DEVICES SPECTRAMAX ID3 MULTI-MODE MICROPLATE READER PROMOTION

Special Introductory Offer

SpectraMax iD3 Multi-Mode Microplate Reader

The SpectraMax® iD3 Multi-Mode Microplate Reader is the cornerstone of a complete laboratory solution to help you expand the boundaries of your research capabilities.

This feature rich reader has a large, high resolution touchscreen with an embedded software package allowing you to quickly set up and run custom protocols; integrates near field communication (NFC) functionality for instrument personalisation; highlights an ultra cooled photomultiplier tube (PMT) that reduces background noise; and has SmartInject™ Technology ensuring equal injector mixing across the plate for high precision experiments. With optimised reagents and the industry leading data acquisition and analysis tool SoftMax® Pro 7 Software, the SpectraMax iD3 allows you to customise your workflow to perfectly match your needs.

MOLECULAR DEVICES

NEW!



Ref. Code: ID3-STD

Valid until 21/12/2017

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
SpectraMax iD3 Plate Reader	iD3 - STD	\$72,920.00	\$49,650.00

SOURCE

PANASONIC DUAL HEAT STERILISATION INCUBATORS

NEW!

New Dual Heat Sterilisation Incubators from Panasonic

Dual Heat Sterilisation

Dual heat sterilisation utilises the incubator's two heaters during the 180°C sterilisation process, which takes 11 hours. Because there is no effect on temperature inside stacked incubators due to low heat dissipation, cell culturing can continue uninterrupted.

Precise & Regulated Environment

InCu-saFe and SafeCell UV both function to prevent contamination. Direct Heat System and melamine foam insulation ensure optimal temperature distribution throughout the chamber while the Dual IR sensor controls the CO₂ level.

Improved Use & Maintenance

A colour LCD touchscreen panel allows full control, even with gloved hands. Transfer of data is easy via a USB port. The easy-to-clean incubator interior features fully rounded corners and integrated shelf supports.



Panasonic

Ref. Code: PANA-NEW

Description	Code	List Price (excl. GST)
CO2 Incubator, 170L, Dual Heat Sterilisation	SANYMCO-170AICD-PE	\$11,499.00
CO2 Incubator, 170L, UV, Dual Heat Sterilisation	SANYMCO-170AICUVD-PE	\$11,949.00

VACUUBRAND, JAVAC & INTEGRA PUMP PROMOTION

UP TO 10% OFF!

Save now on the entire range of vacuum pumps from Vacuubrand and Javac, also included the DOSE IT peristaltic pump from Integra

VACUUBRAND has a 50-year reputation as the innovative expert in vacuum for science. Chemistry diaphragm pumps from VACUUBRAND feature uncompromising chemistry designs. Their construction with fluoropolymers makes them very resistant to chemical vapors from inlet to exhaust and very tolerant to condensates.

For laboratory applications, Javac Vector LD rotary vane vacuum pumps offers high performance with the safety of Automatic System Isolation (ASI) and 3 years' warranty. Vector LT features a Teflon coating on critical parts of the pump, offering greater longevity in aggressive environments.

Integra DOSE IT has been especially designed for laboratories looking for an easy-to-use, compact and portable peristaltic pump. With its light weight and compact design, DOSE IT fits everywhere in the lab and can be easily moved.



Ref. Code: PUMP-PROMO

Valid until 31/12/2017

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
Vacuubrand VARIO PC 3001 Chemistry Pump Unit, speed 2 m ³ /hr, final pressure 2 mbar	VACU696700	\$11,560.00	\$10,404.00
Vacuubrand VARIO PC 3002 Chemistry Pump Unit, speed 2.8 m ³ /hr, final pressure 7 mbar	VACU733500	\$12,300.00	\$11,070.00
Vacuubrand ME1C Chemistry Diaphragm Pump, speed 0,7 m ³ /hr, ultimate vacuum 100 mbar	VACU721100	\$1,850.00	\$1,665.00
Vacuubrand ME2C Chemistry Diaphragm Pump, speed 2.1 m ³ /hr, ultimate vacuum 70 mbar	VACU730100	\$3,030.00	\$2,727.00
Javac Rotary Vane Pump Laboratory 2 stage, 9m ³ /hr, 5 x 10 ⁻⁴ mbar Ultimate Vacuum	JAVAVLD-9-2	\$2,899.00	\$2,610.00
Javac Rotary Vane pump Teflon Laboratory 2 stage, 9m ³ /hr, 5 x 10 ⁻⁴ mbar Ultimate Vacuum	JAVAVLT-9-2	\$3,643.00	\$3,279.00
Integra Peristaltic Pump DOSE IT Complete System, 0,6 ml/min – 5L/min, 0,1ml – 9999ml	INTG171000	\$5,560.00	\$5,004.00

15% OFF SHIMADZU BALANCES



15%
OFF

Save 15% on the entire range of Shimadzu Balances

Shimadzu have been a market leader in balance technology since the early 1900's. More recently, they have introduced the patented UniBloc technology. The UniBloc weighing technology ensures stable temperature characteristics, excellent response time, and stable corner-load performance. In addition, the UniBloc design permits a consistency of production that assures reliability and a long operation life. Balances range from 0.01mg up to 52kg capacities.

Speak to your local Bio-Strategy Equipment Specialist for further details or call Customer Service on 1300 727 696 or email sales.au@bio-strategy.com



Valid until 31/12/2017

Ref. Code: SHIM-SALE

PRIORCLAVE QCS SERIES FRONT LOADING ELECTRICALLY HEATED AUTOCLAVES



Versatile affordable sterilising solutions - from 100 to 500L

Priorclave QCS series Autoclaves with cylindrical chambers and 2 shelves offer versatility and cost efficiency. Constructed to British Standards, the units have grade 316L stainless steel chambers and are inherently safe with thermal interlocks and pressure interlocks. With a choice of 5 capacities – there will be a Priorclave to handle your sterilisation requirements.

- 6 capacities: 100/150/200/320 & 500L
- Single handed lever action door mechanism
- USB port, fan forced cooling and timed free-steaming
- 10 year door and chamber warranty
- Many options and accessories to select from

Speak to your local Bio-Strategy Specialist for further details or call Customer Service on 1800 008 453 or email sales.au@bio-strategy.com



Ref. Code: PRIO-QCS

25% OFF ON ALL HAMILTON PRP HPLC COLUMNS



25%
OFF

Special Introductory Offer

Hamilton offers a complete selection of polymer-based HPLC columns for Reversed-Phase, Anion Exchange, Cation Exchange and Ion Exclusion as well as silica-based columns for Reversed-Phase.

- Better handling of pressure (> 5,000 psi)
- Longer life (can withstand harsher temperature and pH)
- Works better at higher temperatures
- Better sample recovery compared to silica and can be cleaned
- Easily determines anions in almost any sample matrix
- Compatible with any HPLC or ion chromatograph



Ref Code: HAM-HPLC

Valid until 31/12/2017

Description	Code	Sale Price (excl. GST)
PRP-X100 10µm Column 100 x 4.1mm	79439	\$930.00
PRP-X100 10µm Column 150 x 4.1mm	79434	\$1,005.00

SOURCE

HAMILTON ELISA NIMBUS

Designed using the Microlab NIMBUS liquid handling platform, ELISA NIMBUS is a complete walk-away solution for automated ELISA processing

Comprehensive yet compact, the small footprint of ELISA NIMBUS is balanced for sample, reagent and disposable tip capacity, and is capable of processing up to 8 x 96-well microplates in parallel during a run. With 4 independent pipetting channels, a microplate washer, Hamilton Incubator Shaker (HIS), an absorbance reader, and an intuitive, easy-to-use GUI with touchscreen and real-time run control dashboard, the ELISA NIMBUS provides Hamilton quality at a competitive price.

HAMILTON
THE MEASURE OF EXCELLENCE™

ELISA
MADE
EASY!



Ref. Code: HAM-ELISA

Key Features:

- Up to 8 ELISA plates can be run simultaneously
- Sample/Control capacity up to 192 per run
- Complete tracking through barcode reading of sample tubes, sample plates control tubes, reagent troughs and assay plates
- Microplate Reader with preinstalled wavelengths of 405nm, 450nm, 492nm and 620nm as standard and optional 340-750nm also available
- 8-channel Microplate Strip Washer with 3 x 2L wash buffer supply bottles
- One Incubator tower with 4 independent shaking incubators which can be individually temperature controlled from ambient +4°C to 60°C

HAMILTON ADVANCED DILUTER WITH DISPOSABLE TIP HAND PROBE

\$8,775.00 for a limited time only (normally \$11,395.00)

An excellent tool for the applications that require a higher level of assurance.

The sample is aspirated into a disposable plastic tip ("ClickSure") which is thrown away between each sample, eliminating any chance for carryover.

- Aliquot dispensing, Serial dispensing, Dilution, Dry weight dilution, Pipetting and Titration
- Dramatically increases throughput and consistency
- Simple to integrate into a larger process
- Semi-automated and programmable with a touchscreen
- Chemically inert fluid path
- Compliance to EPA, FDA (GLP, GMP) and ISO regulations

Speak to your local Bio-Strategy Specialist for further details or call Customer Service on 1800 008 453 or email sales.au@bio-strategy.com

Ref. Code: HAM-ML625

HAMILTON
THE MEASURE OF EXCELLENCE™

SPECIAL



Valid until 31/12/2017

NEW GE HEALTHCARE ÄKTA CONSUMABLES AND COLUMNS GUIDE

The new ÄKTA consumables and trouble-shooting guide now available in paper or PDF. GE Healthcare products for your ÄKTA chromatography systems & protein sample prep.

Request your ÄKTA Consumables Guide today by calling 1300 727 696 or email sales.au@bio-strategy.com



MP BIO POWERFUL FASTPREP HOMOGENISER WITH LYSING KITS ON SPECIAL!



38% OFF

Grind, Homogenise and Lyse difficult samples in 40 Seconds or Less!
38% discount on FastPrep 5G system and savings on Lysing kits!


Obtain Intact DNA, RNA and proteins from a wide range of sample types such as plant, food, soil, bacteria, animal and human tissue and many others.

Highest speed in combination with its unique multidirectional beating of specialised Lysis Beads on the sample material. The results are highest yields and purity. The FastPrep grinds, homogenises and lyses cell and tissue in 40 seconds or less!

- **Flexibility in sample tube processing sizes:** (24 or 48x2ml, 2x50ml, 12x15ml, 24x4.5ml tubes) at either cryogenic or ambient temperature
- **Exclusive software & touchscreen:** Intuitive, easy touch screen programming or choose from 70 pre-programmed protocols, USB port included
- **High performance:** Increased power and a unique optimised motion lyses samples by multidirectional, simultaneous impact with lysing matrix particles
- **Many years of lysing expertise bundled in one instrument:** Large range of lysing tubes and library of protocols for difficult to lyse samples available! Ask for 10% discount on FastDNA lysing and spin kits for soil, fecal, plant and animal tissue
- **Trade in your old homogeniser and save an additional \$600:** Get the FastPrep 5G for \$11,900 when trading in your old homogeniser (regular List Price: \$19,290, non-trade-in promo price: \$12,500)



Very Simple Workflow



For more details, download or request the brochure

Ref. Code: MPB-FASTPREP

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
FastPrep-24™ 5G Instrument << Trade in price: \$11,900.00 excl. GST >>	1 Unit	ICNA116005500	\$19,290.00	\$12,500.00
Lysing Matrix A Grind or homogenize tough and resistant samples	50x2ml	ICNA116910050	\$316.74	\$285.06
Lysing Matrix A	100x2ml	ICNA116910100	\$447.36	\$402.62
Lysing Matrix A	500x2ml	ICNA116910500	\$993.05	\$893.75
Lysing Matrix A	50x15ml	ICNA116930050	\$452.58	\$407.32
Lysing Matrix A	10x50ml	ICNA116950010	\$489.35	\$440.42
Lysing Matrix B Gram + and - bacteria, fungal tissue and spores	50x2ml	ICNA116911050	\$306.40	\$275.76
Lysing Matrix B	100x2ml	ICNA116911100	\$380.31	\$342.28
Lysing Matrix B	500x2ml	ICNA116911500	\$1,397.53	\$1,257.78
Lysing Matrix B	25x15ml	ICNA116951025	\$221.22	\$199.10
Lysing Matrix B	10x50ml	ICNA116951010	\$225.58	\$203.02
Lysing Matrix C Lysing yeast, fungi and spores, mould and algae	100x2ml	ICNA116912100	\$300.91	\$270.82
Lysing Matrix C	25x15ml	ICNA116932025	\$242.46	\$218.21
Lysing Matrix D Softer tissues, brain, liver, kidney, lung, plants	100x2ml	ICNA116913100	\$395.87	\$356.28
Lysing Matrix D	500x2ml	ICNA116913500	\$723.58	\$651.22
Lysing Matrix D	50x15ml	ICNA116933050	\$584.69	\$526.22
Lysing Matrix E Environmental samples: Soil, sludge, wastewater and feces	50x2ml	ICNA116914050	\$273.81	\$246.43
Lysing Matrix E Mixed Samples: Tissue and RNA extraction	100x2ml	ICNA116914100	\$507.48	\$456.73
Lysing Matrix E	500x2ml	ICNA116914500	\$1,615.22	\$1,453.70

SOURCE

PALL CENTRIFUGAL DEVICES

15% off the Pall Centrifugal Devices range

These devices are used for simple, reliable concentrating and desalting of samples.

Select the low protein binding Omega™ membrane for ultrafiltration applications. The Bio-Inert® membrane is ideal for microfiltration of biological solutions. The GHP membrane is ideal for sample preparation in analytical chemistry (HPLC, GC, IC and MS). It is suitable for aqueous, organic and aggressive organic samples and for rapid processing of samples.

- Rapid processing of samples
- Typical recoveries are greater than 90%
- Housing: Low binding PP
- Ultrasonically welded seals prevent bypass or seal failure



15%
OFF



Ref. Code: PALL-CENT

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
Nanosep®				
Nanosep® MF with Bio-Inert® membrane, 0,2µm, Turquoise	1 * 24 items	PALLODM02C33	\$83.00	\$70.55
Nanosep® MF with Bio-Inert® membrane, 0,2µm, Turquoise	1 * 100 items	PALLODM02C34	\$290.00	\$246.50
Nanosep® MF with Bio-Inert® membrane, 0,2µm, Turquoise	1 * 500 items	PALLODM02C35	\$1,444.00	\$1,227.40
Nanosep® MF with Bio-Inert® membrane, 0,45µm, Purple	1 * 24 items	PALLODM45C33	\$83.00	\$70.55
Nanosep® MF with Bio-Inert® membrane, 0,45µm, Purple	1 * 100 items	PALLODM45C34	\$290.00	\$246.50
Nanosep® MF with Bio-Inert® membrane, 0,45µm, Purple	1 * 500 items	PALLODM45C35	\$1,242.00	\$1,055.70
Nanosep® MF with GHP membrane, 0,45µm, Clear	1 * 100 items	PALLODGHPC34	\$401.00	\$340.85
Nanosep® MF with GHP membrane, 0,45µm, Clear	1 * 500 items	PALLODGHPC35	\$1,508.00	\$1,281.80
Nanosep® with Omega™ membrane, 10kD, Blue	1 * 24 items	PALLOD010C33	\$117.00	\$99.45
Nanosep® with Omega™ membrane, 10kD, Blue	1 * 100 items	PALLOD010C34	\$508.00	\$431.80
Nanosep® with Omega™ membrane, 10kD, Blue	1 * 500 items	PALLOD010C35	\$1,590.00	\$1,351.50
Nanosep® with Omega™ membrane, 100kD, Clear	1 * 24 items	PALLOD100C33	\$113.00	\$96.05
Nanosep® with Omega™ membrane, 100kD, Clear	1 * 100 items	PALLOD100C34	\$379.00	\$322.15
Nanosep® with Omega™ membrane, 100kD, Clear	1 * 500 items	PALLOD100C35	\$1,587.00	\$1,348.95
Nanosep® with Omega™ membrane, 3kD, Grey	1 * 24 items	PALLOD003C33	\$112.00	\$95.20
Nanosep® with Omega™ membrane, 3kD, Grey	1 * 100 items	PALLOD003C34	\$379.00	\$322.15
Nanosep® with Omega™ membrane, 3kD, Grey	1 * 500 items	PALLOD003C35	\$1,587.00	\$1,348.95
Nanosep® with Omega™ membrane, 30kD, Red	1 * 24 items	PALLOD030C33	\$112.00	\$95.20
Nanosep® with Omega™ membrane, 30kD, Red	1 * 100 items	PALLOD030C34	\$379.00	\$322.15
Nanosep® with Omega™ membrane, 30kD, Red	1 * 500 items	PALLOD030C35	\$1,587.00	\$1,348.95
Nanosep® with Omega™ membrane, 300kD, Orange	1 * 24 items	PALLOD300C33	\$113.00	\$96.05
Nanosep® with Omega™ membrane, 300kD, Orange	1 * 100 items	PALLOD300C34	\$379.00	\$322.15
Nanosep® with Omega™ membrane, 300kD, Orange	1 * 500 items	PALLOD300C35	\$1,587.00	\$1,348.95



PALL CENTRIFUGAL DEVICES (CONT.)


**15%
OFF**

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
Microsep™				
Microsep™ Advance with Omega™ membrane, 1kD, Yellow	1 * 24 items	PALLMCP001C41	\$150.00	\$127.50
Microsep™ Advance with Omega™ membrane, 1kD, Yellow	1 * 100 items	PALLMCP001C46	\$540.00	\$459.00
Microsep™ Advance with Omega™ membrane, 10kD, Blue	1 * 24 items	PALLMCP010C41	\$148.00	\$125.80
Microsep™ Advance with Omega™ membrane, 10kD, Blue	1 * 100 items	PALLMCP010C46	\$534.00	\$453.90
Microsep™ Advance with Omega™ membrane, 100kD, Clear	1 * 24 items	PALLMCP100C41	\$148.00	\$125.80
Microsep™ Advance with Omega™ membrane, 100kD, Clear	1 * 100 items	PALLMCP100C46	\$534.00	\$453.90
Microsep™ Advance with Omega™ membrane, 3kD, Grey	1 * 24 items	PALLMCP003C41	\$148.00	\$125.80
Microsep™ Advance with Omega™ membrane, 3kD, Grey	1 * 100 items	PALLMCP003C46	\$534.00	\$453.90
Microsep™ Advance with Omega™ membrane, 30kD, Red	1 * 24 items	PALLMCP030C41	\$148.00	\$125.80
Microsep™ Advance with Omega™ membrane, 30kD, Red	1 * 100 items	PALLMCP030C46	\$534.00	\$453.90
Microsep™ Advance with Supor® membrane, 0,2µm, Turquoise	1 * 24 items	PALLMCPM02C67	\$148.00	\$125.80
Microsep™ Advance with Supor® membrane, 0,2µm, Turquoise	1 * 100 items	PALLMCPM02C68	\$540.00	\$459.00
Microsep™ Advance with Supor® membrane, 0,45µm, Purple	1 * 24 items	PALLMCPM45C67	\$148.00	\$125.80
Microsep™ Advance with Supor® membrane, 0,45µm, Purple	1 * 100 items	PALLMCPM45C68	\$534.00	\$453.90
Macrosep®				
Macrosep® Advance with Omega™ membrane, 1kD, Yellow	1 * 6 items	PALLMAP001C36	\$78.00	\$66.30
Macrosep® Advance with Omega™ membrane, 1kD, Yellow	1 * 24 items	PALLMAP001C37	\$279.00	\$237.15
Macrosep® Advance with Omega™ membrane, 1kD, Yellow	1 * 100 items	PALLMAP001C38	\$1,133.00	\$963.05
Macrosep® Advance with Omega™ membrane, 10kD, Blue	1 * 6 items	PALLMAP010C36	\$77.00	\$65.45
Macrosep® Advance with Omega™ membrane, 10kD, Blue	1 * 24 items	PALLMAP010C37	\$276.00	\$234.60
Macrosep® Advance with Omega™ membrane, 10kD, Blue	1 * 100 items	PALLMAP010C38	\$1,081.00	\$918.85
Macrosep® Advance with Omega™ membrane, 100kD, Clear	1 * 6 items	PALLMAP100C36	\$77.00	\$65.45
Macrosep® Advance with Omega™ membrane, 100kD, Clear	1 * 24 items	PALLMAP100C37	\$276.00	\$234.60
Macrosep® Advance with Omega™ membrane, 100kD, Clear	1 * 100 items	PALLMAP100C38	\$1,081.00	\$918.85
Macrosep® Advance with Omega™ membrane, 3kD, Grey	1 * 6 items	PALLMAP003C36	\$77.00	\$65.45
Macrosep® Advance with Omega™ membrane, 3kD, Grey	1 * 24 items	PALLMAP003C37	\$276.00	\$234.60
Macrosep® Advance with Omega™ membrane, 3kD, Grey	1 * 100 items	PALLMAP003C38	\$1,081.00	\$918.85
Macrosep® Advance with Omega™ membrane, 30kD, Red	1 * 6 items	PALLMAP030C36	\$77.00	\$65.45
Macrosep® Advance with Omega™ membrane, 30kD, Red	1 * 24 items	PALLMAP030C37	\$276.00	\$234.60
Macrosep® Advance with Omega™ membrane, 30kD, Red	1 * 100 items	PALLMAP030C38	\$1,081.00	\$918.85
Macrosep® Advance with Supor® membrane, 0,2µm, Turquoise	1 * 24 items	PALLMAPM02C67	\$276.00	\$234.60
Macrosep® Advance with Supor® membrane, 0,2µm, Turquoise	1 * 100 items	PALLMAPM02C68	\$1,081.00	\$918.85
Macrosep® Advance with Supor® membrane, 0,45µm, Purple	1 * 24 items	PALLMAPM45C67	\$276.00	\$234.60
Macrosep® Advance with Supor® membrane, 0,45µm, Purple	1 * 100 items	PALLMAPM45C68	\$1,081.00	\$918.85
Jumbosep™				
Jumbosep™ Membrane insert, 10kD, Blue	1 * 12 items	PALLOD010C65	\$289.00	\$245.65
Jumbosep™ Membrane insert, 100kD, Clear	1 * 12 items	PALLOD100C65	\$291.00	\$247.35
Jumbosep™ Membrane insert, 3kD, Grey	1 * 12 items	PALLOD003C65	\$289.00	\$245.65
Jumbosep™ Membrane insert, 30kD, Red	1 * 12 items	PALLOD030C65	\$289.00	\$245.65
Jumbosep™ Membrane insert, 300kD, Orange	1 * 12 items	PALLOD300C65	\$305.00	\$259.25



SOURCE

30% OFF AZLON MEASURING CYLINDERS

Class B Measuring Cylinders on sale!

Polypropylene Class B Measuring Cylinders with large pouring spout for easy dispensing.

- Stable hexagonal base with anti-suction design that prevents base from sticking to a wet surface
- Permanent moulded graduations
- Good chemical resistance and extremely durable
- Recycling codes and batch code on each cylinder
- Autoclavable
- ISO 6706



30%
OFF

Azlon

Ref. Code: AZLO-CYLIN

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
Measuring Cylinder, Spouted with moulded Graduations, 25ml	1 * 10 items	AZLOCP0025P	\$36.90	\$25.83
Measuring Cylinder, Spouted with moulded Graduations, 50ml	1 * 10 items	AZLOCP0050P	\$43.70	\$30.59
Measuring Cylinder, Spouted with moulded Graduations, 100ml	1 * 5 items	AZLOCP0100P	\$28.70	\$20.09
Measuring Cylinder, Spouted with moulded Graduations, 250ml	1 * 1 items	AZLOCP0250P	\$8.70	\$6.09
Measuring Cylinder, Spouted with moulded Graduations, 500ml	1 * 1 items	AZLOCP0500P	\$12.60	\$8.82
Measuring Cylinder, Spouted with moulded Graduations, 1000ml	1 * 1 items	AZLOCP1000P	\$21.10	\$14.77
Measuring Cylinder, Spouted with moulded Graduations, 2000ml	1 * 1 items	AZLOCP2000P	\$42.80	\$29.96

40% OFF PYREX PREMIUM VOLUMETRIC FLASKS!

PYREX® Volumetric Flasks, borosilicate glass with polyethylene stopper

- Calibrated to the tolerances required by the United States Pharmacopeia (USP), which are tighter than those of ISO/DIN standards
- Each flask is laser etched with a unique individual serial number, batch code and date of manufacture - Batch certificates are downloadable from website
- Conformity 'H', certified
- Manufactured from chemically resistant Pyrex® glass
- White enamel inscriptions
- BS EN ISO 1042, DIN 12664



40%
OFF


SciLabware
AZLON® MBL® PYREX® QUICKFIT® WHEATON®

Ref. Code: BIBB-FLASK

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
Pyrex Volumetric Flask, Class A with stopper, 1ml	1 * 5 items	BIBB3151/01	\$323.00	\$193.80
Pyrex Volumetric Flask, Class A with stopper, 5ml	1 * 5 items	BIBB3151/04	\$81.20	\$48.72
Pyrex Volumetric Flask, Class A with stopper, 10ml	1 * 5 items	BIBB3151/06	\$82.90	\$49.74
Pyrex Volumetric Flask, Class A with stopper, 20ml	1 * 5 items	BIBB3151/08	\$85.30	\$51.18
Pyrex Volumetric Flask, Class A with stopper, 50ml	1 * 5 items	BIBB3151/12	\$95.80	\$57.48
Pyrex Volumetric Flask, Class A with stopper, 100ml	1 * 5 items	BIBB3151/14	\$96.20	\$57.72
Pyrex Volumetric Flask, Class A with stopper, 200ml	1 * 5 items	BIBB3151/16	\$146.00	\$87.60
Pyrex Volumetric Flask, Class A with stopper, 250ml	1 * 2 items	BIBB3151/18	\$66.70	\$40.02
Pyrex Volumetric Flask, Class A with stopper, 500ml	1 * 2 items	BIBB3151/20	\$77.20	\$46.32
Pyrex Volumetric Flask, Class A with stopper, 1000ml	1 * 2 items	BIBB3151/22	\$111.00	\$66.60
Pyrex Volumetric Flask, Class A with stopper, 2000ml	1 * 1 items	BIBB3151/24	\$111.00	\$66.60
Pyrex Volumetric Flask, Class A with stopper, 5000ml	1 * 1 items	BIBB3151/26	\$460.00	\$276.00

LABADVANTAGE MICROCENTRIFUGE TUBES

SAVE
20%

Give your Lab the Advantage, the LabAdvantage!

Made from polypropylene, these tubes are graduated and autoclavable

- Superior clarity and strength
- Temperature rating from -80°C to 121°C
- Fit all standard rotors
- Positive seal design allows for repeated cap closures and eliminates pressure induced popping of caps
- Can withstand up to 26,000xg
- Assorted colour pack includes one each of blue, green, yellow, red and purple tubes



Ref. Code: SORE-MCT

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
Microcentrifuge tube with attached lid, 0.65ml, Natural	1 * 1,000 items	SORELAC11154	\$32.70	\$26.16
Microcentrifuge tube with attached lid, 0.65ml, Assorted colours	1 * 1,000 items	SORELAC11214	\$35.40	\$28.32
Microcentrifuge tube with attached lid, low binding, 0.65ml, Colourless	1 * 500 items	SORELAC1130	\$32.70	\$26.16
Microcentrifuge tube with attached lid, 1.7ml, Blue	1 * 500 items	SORELAC11434	\$18.00	\$14.40
Microcentrifuge tube with attached lid, 1.7ml, Colourless	1 * 500 items	SORELAC11514	\$16.50	\$13.20
Microcentrifuge tube with attached lid, 1.7ml, Assorted colours	1 * 500 items	SORELAC11594	\$19.10	\$15.28
Microcentrifuge tube with attached lid, 2.0ml, Colourless	1 * 400 items	SORELAC12004	\$16.50	\$13.20
Microcentrifuge tube with attached lid, 2.0ml, Assorted colours	1 * 400 items	SORELAC12064	\$19.10	\$15.28
Microcentrifuge tube with attached lid, Dolphin nose tip, 2.0ml, Natural	1 * 5,000 items	SORELAC2153	\$261.00	\$208.80
Microcentrifuge tube with attached lid, low binding, 1.7ml, Colourless	1 * 250 items	SORELAC3964T	\$17.00	\$13.60

25% OFF KARTELL CUVETTES

Specials on Spectro Cuvettes

Cuvettes for spectrophotometry application in the visible and UV spectrum. Clear and sharp optical path with indication of the optical path direction. They are made of optical PS and UV grade PMMA and are equipped with optically machined walls (optical faces) to ensure optimal transparency along the spectral field from 340 to 800nm (optical PS) and from 280 to 800nm (UV grade PMMA).

Available in the macro version (4.5ml), low-form semi-micro (1.5ml), high-form semi-micro (2.5ml), with 2 or 4 optical faces. The two-optical face version has grooves on the two walls not optically machined to provide a better grip and to prevent positioning errors.

Packaged in dust-free conditions.

Kartell

25%
OFF!



Ref. Code: KART-CUVE

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
Cuvette, UV, 4 clear sides PMMA, 4.50ml, 10mm	1 * 100 items	KART1961	\$23.10	\$17.33
Cuvette, UV, Macro PMMA, 4.50ml, 10mm	1 * 100 items	KART1939	\$26.80	\$20.10
Cuvette, UV, Semi-micro PMMA, 1.50ml, 10mm	1 * 100 items	KART1948	\$50.30	\$37.73
Cuvette, UV, Semi-micro PMMA, 2.50ml, 10mm	1 * 100 items	KART1941	\$51.20	\$38.40
Cuvette, VIS, 4 clear faces PS, 4.50ml, 10mm	1 * 100 items	KART1960	\$30.90	\$23.18
Cuvette, VIS, Macro PS, 4.50ml, 10mm	1 * 100 items	KART1937	\$37.00	\$27.75
Cuvette, VIS, Semi-micro PS, 1.50ml, 10mm	1 * 100 items	KART1938	\$30.60	\$22.95
Cuvette, VIS, Semi-micro PS, 2.50ml, 10mm	1 * 100 items	KART1940	\$25.80	\$19.35

SOURCE

15% OFF VWR TRANSFER PIPETTES

Great saving on these VWR Transfer Pipettes

VWR Transfer Pipettes made from LDPE, disposable

These unbreakable all-in-one pipettes eliminate the potential hazard of broken glass and exposure to infectious materials. The low-affinity surface reduces the loss of cells and valuable proteins due to binding. They work well whenever there is a need for quick, safe transfer of fluids.

- Will not shatter or break
- Non-toxic and inert to biological fluids and most acids
- Can be sealed and refrigerated; sterilised and non-sterilised available
- No bulb to insert or remove
- Uniform drop size



15%
OFF



Ref. Code: VWR-PIPE

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
23.0ml, Non graduated, 300mm, Non-Sterile	1 * 100 items	VWRI612-3627	\$7.50	\$6.53
7.0ml, Graduated, large bulb, 153mm, Non-Sterile	1 * 500 items	VWRI612-4472	\$23.10	\$20.10
5.8ml, Non graduated, fine tip, standard, 147mm, Non-Sterile	1 * 500 items	VWRI612-4492	\$20.90	\$18.18
5.8ml, Graduated, large bulb, 155mm, Non-Sterile	1 * 500 items	VWRI612-4494	\$22.10	\$19.23
5.8ml, Graduated, large bulb, 155mm, Sterile	1 * 500 items	VWRI612-4497	\$23.40	\$20.36
7.5ml, Non graduated, standard, 155mm, Non-Sterile	1 * 500 items	VWRI612-4504	\$20.60	\$17.92
3.0ml, Graduated, small bulb, 140mm, Non-Sterile	1 * 500 items	VWRI612-4520	\$16.30	\$14.18
7.5ml, Graduated, large bulb, 155mm, Non-Sterile	1 * 500 items	VWRI612-4545	\$22.00	\$19.14
1.2ml, Non graduated, narrow stem, short, 64mm, Non-Sterile	1 * 500 items	VWVU16001-170	\$34.00	\$29.58
1.7ml, Non graduated, general purpose, 87mm, Non-Sterile	1 * 500 items	VWVU16001-174	\$20.90	\$18.18
3.0ml, Graduated, small bulb, 138mm, Non-Sterile	1 * 500 items	VWVU16001-176	\$14.60	\$12.70
5.0ml, Graduated, large bulb, 155mm, Non-Sterile	1 * 500 items	VWVU16001-180	\$14.60	\$12.70
5.0ml, Graduated, large bulb, 155mm, Non-Sterile	1 * 500 items	VWVU16001-182	\$14.60	\$12.70
7.5ml, Graduated, large bulb, 155mm, Non-Sterile	1 * 500 items	VWVU16001-188	\$14.60	\$12.70
1.5ml, Non graduated, extended fine tip, small bulb, 104mm, Non-Sterile	1 * 400 items	VWVU16001-192	\$14.60	\$12.70
23.0ml, Non graduated, 300mm, Non-Sterile	1 * 100 items	VWVU16001-196	\$6.70	\$5.83
5.0ml, Non graduated, general purpose, 155mm, Non-Sterile	1 * 500 items	VWVU414004-001	\$14.60	\$12.70
7.5ml, Non graduated, standard, 155mm, Sterile	1 * 500 items	VWVU414004-005	\$32.80	\$28.54
5.8ml, Graduated, large bulb, 155mm, Non-Sterile	1 * 500 items	VWVU414004-014	\$14.60	\$12.70
5.8ml, Graduated, large bulb, 155mm, Sterile	1 * 500 items	VWVU414004-015	\$19.70	\$17.14
3.0ml, Graduated, small bulb, 140mm, Non-Sterile	1 * 500 items	VWVU414004-017	\$14.60	\$12.70
5.8ml, Non graduated, fine tip, standard, 147mm, Non-Sterile	1 * 500 items	VWVU414004-020	\$14.60	\$12.70
5.8ml, Non graduated, fine tip, standard, 147mm, Sterile	1 * 500 items	VWVU414004-021	\$94.50	\$82.22
5.0ml, Non graduated, extended fine tip, large bulb, 152mm, Non-Sterile	1 * 250 items	VWVU414004-023	\$10.70	\$9.31
1.5ml, Graduated, thin stem, 116mm, Sterile	1 * 500 items	VWVU414004-024	\$32.80	\$28.54
1.5ml, Graduated, thin stem, 116mm, Sterile	1 * 500 items	VWVU414004-025	\$15.90	\$13.83
5.5ml, Non graduated, 229mm, Sterile	1 * 400 items	VWVU414004-026	\$123.00	\$107.01
5.5ml, Non graduated, 229mm, Sterile	1 * 400 items	VWVU414004-028	\$17.20	\$14.96
9.3ml, Non graduated, 138mm, Non-Sterile	1 * 250 items	VWVU414004-032	\$8.10	\$7.05
5.8ml, Graduated, large bulb, 155mm, Sterile	1 * 500 items	VWVU414004-040	\$15.90	\$13.83
7.5ml, Graduated, large bulb, 155mm, Sterile	1 * 500 items	VWVU414004-045	\$17.20	\$14.96
5.8ml, Non graduated, fine tip, standard, 147mm, Non-Sterile	1 * 500 items	VWVU414004-048	\$19.70	\$17.14

\$7 BOLLÉ AXIS 2 SAFETY SPECTACLES

Safety without compromising comfort!

Bollé AXIS 2 offers protection adapted to suit each user thanks to its pivoting temple arms

- High contour (10.5 base) single piece lens
- Co-moulded tilt adjustable temple arms for personalised and secure fit and protection against dust and other particles
- Adjustable silicone nose pads for comfortable fit
- Fit: Unisex - Medium / Large
- Anti-scratch/Anti-fog coating
- Approved to AS/NZS 1337.1.2010 for Medium Impact Protection

bollé
SAFETY

Only
\$7!



Ref. Code: **BOLE-AXIS2**

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
AXIS 2 Safety Spectacle, adjustable arms and nose bridge, clear lens	1 * 1 items	BOLE1654101	\$8.00	\$7.00
AXIS 2 Safety Spectacle, adjustable arms and nose bridge, smoke lens	1 * 1 items	BOLE1654102	\$8.00	\$7.00
AXIS 2 Safety Spectacle, adjustable arms and nose bridge, contrast lens	1 * 1 items	BOLE1654118	\$8.00	\$7.00

12% OFF ANSELL AFFINITY NEOPRENE GLOVES

Buy now and save!

Micro-Touch® Affinity® Neoprene, powder-free glove with textured fingertips

These examination gloves with beaded cuffs deliver the protection, comfort, stretchability, fit, and tactile sensitivity of natural rubber latex, but are powder and latex-free, making them an ideal choice for latex allergy sufferers. Tested for use with chemotherapy drugs.

- Powder-free to eliminate powder related complications
- Ambidextrous gloves with textured surface finish for better grip in wet and dry conditions
- Thiurams-free and mercaptobenzothiazoles-free, to reduce the risk of contact dermatitis
- Green coloured neoprene for easy identification as a synthetic glove
- Easy donning due to polymer coating
- AQL 1.5



Ansell



12%
OFF

Ref. Code: **ANSE-AFFIN**

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
Micro-Touch® Affinity® Neoprene glove - Extra Small	1 * 100 items	ANSE3770	\$17.05	\$15.00
Micro-Touch® Affinity® Neoprene glove - Small	1 * 100 items	ANSE3771	\$17.05	\$15.00
Micro-Touch® Affinity® Neoprene glove - Medium	1 * 100 items	ANSE3772	\$17.05	\$15.00
Micro-Touch® Affinity® Neoprene glove - Large	1 * 100 items	ANSE3773	\$17.05	\$15.00
Micro-Touch® Affinity® Neoprene glove - Extra Large	1 * 100 items	ANSE3774	\$17.05	\$15.00

WINE/BREWERY KITS & REAGENTS FROM MEGAZYME

Setting new standards for purity, quality and innovation in test technology. Since 1989, Megazyme has become a globally recognised leader in the development and supply of innovative diagnostic test kits and analytical reagents for the food, feed, fermentation, dairy and wine industries.

For more information phone 1300 727 696 or email sales.au@bio-strategy.com



SOURCE

BIO-STRATEGY - YOUR ONE STOP SHOP FOR LABORATORY CHEMICALS!

A large choice of High Quality Chemical Brands – All in one order!

Buy \$600 worth of mixed Laboratory Chemicals including VWR Chemicals and get \$100 off!

Chemicals order consolidation means savings:

Bio-Strategy supports laboratories all over Australia and New Zealand with a vast range of high quality chemicals. Our supplier partners feature unique specialities and purity grades that match your research or process steps. We leave the choice to you to fill the stock room with only the best. Experienced laboratory managers and scientists save money by consolidating their orders.

\$100 savings on your next Chemical order:

October to December only: Contact our Customers Service team and mention the Promo Reference Code CHEMS-SALE with a mixed order of \$600 or more of Laboratory Chemicals to get \$100 immediately!*


*All brands shown above are eligible but all orders must include 1/3 VWR Chemicals. Cannot be combined with other discounts or specials.

Top Runners:

We can help you even more with special prices for your top used chemicals. Contact us today with a list of chemicals and ask for YOUR special price!



SAVE \$100!



Ref. Code: CHEMS-SALE

Valid until 31/12/2017

15% OFF VWR BOXES FOR BROKEN GLASS

Bargains on VWR waste boxes!

Corrugated cardboard with polyethylene pouch, 50 µm thick

Ideal for collecting, transporting and removing broken glass from the laboratory. Consists of a sturdy box made of strong corrugated cardboard and a polyethylene insertion pouch which prevents any possible escape of fluids adhering to the broken glass.

- Folded flat for space saving storage
- A flap is built into the lid for easy filling
- The internal safety lid is closed for the removal of the container
- Delivered unassembled, assembly takes a few seconds



15% OFF!



Ref. Code: VWR-BOX

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
Disposal Box Glass/Liquid Floor Standing, 305x305x685mm	1 * 6 items	VWRI119-0036	\$72.00	\$61.20
Disposal Box Glass/Liquid Bench Standing, 203x203x254mm	1 * 6 items	VWRI119-0037	\$50.60	\$43.01

COMING SOON - New combined Bio-Strategy Brand Index!

We are very excited about combining the former VWR Australia and New Zealand businesses with Bio-Strategy and we are looking forward to offering you an increased number of quality brands, application specialists and extended world class services.

For more information or to request a quotation phone 1300 727 696 or email sales.au@bio-strategy.com



GE WHATMAN FOR CHEMISTRY BROCHURE!

Save up to 47% on products from the Whatman Chemistry Laboratory brochure

Request or download the brochure to see a wide range of products ideal for the chemistry laboratory.

For a limited time only, receive up to 47% off products included in the brochure.


Products in the brochure include:

- Benchkote
- Phase Separator Paper
- Syringe Filters
- Indicator Paper
- And much more....



UP TO
47% OFF

Distributor
GE Healthcare



Ref Code: WHAT-CHEM

Valid until 31/12/2017

Description	UoM	Code	Sale Price (excl. GST)
Benchkote Surface Protector, 460x570 mm sheet	1 * 50 items	WHAT2300-916	\$75.00
Cellulose Soxhlet Extraction Thimble, 1 x 22x80 mm	1 * 25 items	WHAT2800-228	\$73.80
pH strips, colour bonded 0.0 to 14.0 range, 6x80 mm	1 * 100 items	WHAT2613-991	\$19.70
Kjeldahl Analysis Weighing Boat, Grade 609, 55x10x10 mm	1 * 100 items	WHAT10313032	\$118.00
Phase separating paper, 1PS for Solvent Extraction, 125 mm circle	1 * 100 items	WHAT2200-125	\$39.40
Syringe filters, GD/F, 25mm, 0,7µm, 4.6 cm ²	1 * 150 items	WHAT6890-2507	\$386.10

More products on sale - download or request the Whatman Chemistry flyer today!


NEW VWR CHEMICALS CATALOGUE 2017-2019

Bio-Strategy - your partner for all needs in research and production chemicals

- Thousands of products for analysis, research and production, many new to the VWR range
- Multiple grades ensuring the right quality for each application
- Laboratory and production pack sizes
- Custom manufacturing services



Make sure your safety officer gets a copy
by requesting one today from your local
Bio-Strategy Sales Representative



SOURCE

CORNING MATRIGEL AND COOLCELL STARTER KITS

Purchase anything from the Corning Matrigel range and receive 15% off the CoolCell Starter Kits!

The Corning® XT Starter is an open-platform ice-free cooler that accommodates most CoolRack® and CoolSink® modules. Low profile and a small footprint make it ideal for use in the hood, keeping samples cold (0,5 to 4,0°C) up to four hours. 1° to 8°C temperature indicator provides visual assurance of temperature performance. To extend the cooling duration, keep an additional XT cooling core in the freezer and rotate the cores as needed prevent positioning errors.



CORNING

15% OFF!



Bundle up and save!

Ref. Code: CORN-COOL

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
Corning® XT Starter Complete, Low profile, Open Platform Ice-free Cooler	1 * 1 items	CORN432014	\$427.00	\$362.95
Corning® XT Starter Holder Only, Open Platform Ice-free Cooler	1 * 1 items	CORN432015	\$133.00	\$113.05
Corning® XT Starter PCR96 Kit for 96-well PCR plate, strip-wells or 200µl tubes	1 * 1 items	CORN432016	\$1,080.00	\$918.00
Corning® XT Starter M-PCR Kit for 12 x 1.5mL or 2.0mL microfuge tubes and 6 rows of PCR stripwells	1 * 1 items	CORN432017	\$1,060.00	\$901.00
CoolSink XT96F, Flat Bottom Plate Module	1 * 1 items	CORN432070	\$608.00	\$516.80

NEW PALL LABORATORY MANIFOLD

A unique microbiology manifold that allows you to optimise testing without sacrificing cleanliness. With interchangeable components and coupling devices, the manifold takes direction from you.

Applications include:


- Bioburden testing
- Detecting spoilage micro-organisms
- Monitoring for compliance of drinking water
- Membrane Filtration Technique testing

For more information contact your local Bio-Strategy Sales Representative or call Customer Service on 1800 727 696 or email sales.au@bio-strategy.com



Contact us for more information and to place your orders

Bio-Strategy Pty Ltd.
T: 1300 727 696
E: sales.au@bio-strategy.com
www.bio-strategy.com



Prices & discounts in this brochure may be subject to change. Promotions are not valid with any other special offer or discounts and cannot be applied in retrospect. All offers subject to availability and while stocks last. All prices are in AUD excluding GST and local shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.