

Eppendorf 5424 and 5424R Centrifuge Promotion

VALID UNTIL
31 December 2017
Ref. No. EPPE-5424

Enquire
today for
other 5424
package
options!

Equipped for all modern Life Science applications - SAVE NOW!

Eppendorf 5424 and 5424 R (refrigerated) centrifuges are the new laboratory standard. With their 24-place capacity and speed up to $21,130 \times g$ they are perfectly equipped for all modern molecular biology applications in Eppendorf and PCR tubes.

- 5424 R centrifuge cools down to 4 °C in only 8 minutes
- 5424 R centrifuge temperature range from -10 °C to +40 °C
- Soft-touch one-finger lid closure for ergonomic operation
- OptiBowl® design for quiet operation, even without a rotor lid
- SOFT brake to protect sensitive samples
- Certified aerosol-tight rotor comes standard for safe centrifugation of hazardous samples

Description	Product Code	Sales Price (ex. GST)	Sales Price (inc GST)
Centrifuge 5424R refrigerated knob version with FA-45-24-11 rotor	EPPE5404000.480	\$5,945.00	\$6,539.50
Centrifuge 5424 knob version with FA-45-24-11 rotor	EPPE5424000.487	\$2,819.00	\$3,100.90



BIO-STRATEGY PTY LTD.

T: 1300 727 696
W: bio-strategy.com
E: sales.au@bio-strategy.com



bio-strategy
delivering technology

Prices & discounts above may be subject to change. Not valid with any other special offer or discounts and cannot be applied in retrospect. All offers subject to availability and while stocks last. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.