

MP BIO POWERFUL FASTPREP HOMOGENISER WITH LYSING KITS ON SPECIAL!





Grind, Homogenise and Lyse difficult samples in 40 Seconds or Less!
38% discount on FastPrep 5G system and savings on Lysing kits!


Obtain Intact DNA, RNA and proteins from a wide range of sample types such as plant, food, soil, bacteria, animal and human tissue and many others.

Highest speed in combination with its unique multidirectional beating of specialised Lysis Beads on the sample material. The results are highest yields and purity. The FastPrep grinds, homogenises and lyses cell and tissue in 40 seconds or less!

- **Flexibility in sample tube processing sizes:** (24 or 48x2ml, 2x50ml, 12x15ml, 24x4.5ml tubes) at either cryogenic or ambient temperature
- **Exclusive software & touchscreen:** Intuitive, easy touch screen programming or choose from 70 pre-programmed protocols, USB port included
- **High performance:** Increased power and a unique optimised motion lyses samples by multidirectional, simultaneous impaction with lysing matrix particles
- **Many years of lysing expertise bundled in one instrument:** Large range of lysing tubes and library of protocols for difficult to lyse samples available! Ask for 10% discount on FastDNA lysing and spin kits for soil, fecal, plant and animal tissue
- **Trade in your old homogeniser and save an additional \$600:** Get the FastPrep 5G for \$11,900 when trading in your old homogeniser (regular List Price: \$19,290, non-trade-in promo price: \$12,500)



Very Simple Workflow



For more details, download or request the brochure

Ref. Code: MPB-FASTPREP

Valid until 31/12/2017

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
FastPrep-24™ 5G Instrument	<< Trade in price: \$11,900.00 excl. GST >>	ICNA116005500	\$19,290.00	\$12,500.00
Lysing Matrix A Grind or homogenize tough and resistant samples	50x2ml	ICNA116910050	\$316.74	\$285.06
Lysing Matrix A	100x2ml	ICNA116910100	\$447.36	\$402.62
Lysing Matrix A	500x2ml	ICNA116910500	\$993.05	\$893.75
Lysing Matrix A	10x15ml	ICNA116930050	\$452.58	\$407.32
Lysing Matrix A	10x50ml	ICNA116950010	\$489.35	\$440.42
Lysing Matrix B Gram + and - bacteria, fungal tissue and spores	50x2ml	ICNA116911050	\$306.40	\$275.76
Lysing Matrix B	100x2ml	ICNA116911100	\$380.31	\$342.28
Lysing Matrix B	500x2ml	ICNA116911500	\$1,397.53	\$1,257.78
Lysing Matrix B	25x15ml	ICNA116951025	\$221.22	\$199.10
Lysing Matrix B	10x50ml	ICNA116951010	\$225.58	\$203.02
Lysing Matrix C Lysing yeast, fungi and spores, mould and algae	100x2ml	ICNA116912100	\$300.91	\$270.82
Lysing Matrix C	25x15ml	ICNA116932025	\$242.46	\$218.21
Lysing Matrix D Softer tissues, brain, liver, kidney, lung, plants	100x2ml	ICNA116913100	\$395.87	\$356.28
Lysing Matrix D	500x2ml	ICNA116913500	\$723.58	\$651.22
Lysing Matrix D	50x15ml	ICNA116933050	\$584.69	\$526.22
Lysing Matrix E Environmental samples: Soil, sludge, wastewater and feces	50x2ml	ICNA116914050	\$273.81	\$246.43
Lysing Matrix E Mixed Samples: Tissue and RNA extraction	100x2ml	ICNA116914100	\$507.48	\$456.73
Lysing Matrix E	500x2ml	ICNA116914500	\$1,615.22	\$1,453.70

BIO-STRATEGY PTY LTD.

T: 1300 727 696
 W: bio-strategy.com
 E: sales.au@bio-strategy.com



Prices & discounts above may be subject to change. Not valid with any other special offer or discounts and cannot be applied in retrospect. All offers subject to availability and while stocks last. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.